

ROTH* Prescription

Maxillary	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Central	12	5	0	3	UR	4881-0111-10	4821-0111-10
					UL	4881-0121-10	4821-0121-10
Lateral	8	9	0	3	UR	4881-0112-10	4821-0112-10
					UL	4881-0122-10	4821-0122-10
Cuspid	0	11	4	3	UR	4881-0113-10	4821-0113-10
					UL	4881-0123-10	4821-0123-10
Cuspid HK	0	11	4	3	UR	4881-0113-11	4821-0113-11
					UL	4881-0123-11	4821-0123-11
1st & 2nd Bicuspid	-7	0	0	3.2	U	4881-0104-10	4821-0104-10
1st & 2nd Bicuspid HK	-7	0	0	3.2	UR	4881-0114-11	4821-0114-11
					UL	4881-0124-11	4821-0124-11
Kit						4881-000K-H*	4821-000K-H*

*indicate Hook on 3,45, 345 when ordering with hooks

ROTH* Prescription

Mandibular	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Anterior	0	0	0	3	U	4881-0201-10	4821-0201-10
Cuspid	-11	5	2	3	LR	4881-0213-10	4821-0213-10
					LL	4881-0223-10	4821-0223-10
Cuspid HK	-11	5	2	3	LR	4881-0213-11	4821-0213-11
					LL	4881-0223-11	4821-0223-11
1st Bicuspid	-17	0	0	3.2	U	4881-0204-10	4821-0204-10
1st bicuspid HK	-17	0	0	3.2	LR	4881-0214-11	4821-0214-11
					LL	4881-0224-11	4821-0224-11
2nd Bicuspid	-22	0	0	3.2	U	4881-0205-10	4821-0205-10
2nd Bicuspid HK	-22	0	0	3.2	LR	4881-0215-11	4821-0215-11
					LL	4881-0225-11	4821-0225-11

MBT* Prescription

Maxillary	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Central	17	4	0	3	UR	4882-0111-20	4822-0111-20
					UL	4882-0121-20	4822-0121-20
Lateral	10	8	0	3	UR	4882-0112-20	4822-0112-20
					UL	4882-0122-20	4822-0122-20
Cuspid	-7	8	0	3	UR	4882-0113-00	4822-0113-00
					UL	4882-0123-00	4822-0123-00
Cuspid HK	0	8	0	3	UR	4882-0113-11	4822-0113-11
					UL	4882-0123-11	4822-0123-11
1st & 2nd Bicuspid	-7	0	0	3.2	U	4882-0104-10	4822-0104-10
1st & 2nd Bicuspid HK	-7	0	0	3.2	UR	4882-0114-11	4822-0114-11
					UL	4882-0124-11	4822-0124-11
Kit						4882-000K-H*	4822-000K-H*

*indicate Hook on 3,45, 345 when ordering with hooks

MBT* Prescription

Mandibular	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Anterior	-6	0	0	3	U	4882-0201-00	4822-0201-00
Cuspid	-6	3	0	3	LR	4882-0213-00	4822-0213-00
					LL	4882-0223-00	4822-0223-00
Cuspid HK	0	3	0	3	LR	4882-0213-11	4822-0213-11
					LL	4882-0223-11	4822-0223-11
1st Bicuspid	-12	0	0	3.2	U	4882-0204-20	4822-0204-20
1st Bicuspid HK	-12	0	0	3.2	LR	4882-0214-21	4822-0214-21
					LL	4882-0224-21	4822-0224-21
2nd Bicuspid	-17	0	0	3.2	U	4882-0205-20	4822-0205-20
2nd Bicuspid	-17	0	0	3.2	LR	4882-0215-21	4822-0215-21
					LL	4882-0225-21	4822-0225-21

AURORA

TWIN CLEAR BRACKET

FEATURES

- Mechanical lock base for increased retention
- Rounded corners for enhanced patient comfort
- Torque-In-Base for level slot lineup and accurate finishing
- Compound Contour Base for an ideal bracket fit

