

ROTH* Prescription

Maxillary	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Central	12	5	0	3.5	UL	4181-0121-10	4121-0121-10
					UR	4181-0111-10	4121-0111-10
Lateral	8	9	0	3	UL	4181-0122-10	4121-0122-10
					UR	4181-0112-10	4121-0112-10
Cuspid	0	11	0	3.3	UL	4181-0123-10	4121-0123-10
					UR	4181-0113-10	4121-0113-10
Cuspid HK	0	11	0	3.3	UL	4181-0123-11	4121-0123-11
					UR	4181-0113-11	4121-0113-11
1st & 2nd Bicuspid	-7	0	2	3	UL	4181-0124-10	4121-0124-10
					UR	4181-0114-10	4121-0114-10
1st & 2nd Bicuspid HK	-7	0	2	3	UL	4181-0124-11	4121-0124-11
					UR	4181-0114-11	4121-0114-11
Kit						4181-000K-H*	4121-000K-H*

*indicate Hook on 3,45, 345 when ordering with hooks

ROTH* Prescription

Mandibular	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Anterior	0	0	0	2.5	L	4181-0201-10	4121-0201-10
Cuspid	-11	5	0	3	LL	4181-0223-10	4121-0223-10
					LR	4181-0213-10	4121-0213-10
Cuspid HK	-11	5	0	3	LR	4181-0213-11	4121-0213-11
					LL	4181-0223-11	4121-0223-11
1st Bicuspid	-17	0	4	3	LL	4181-0224-10	4121-0224-10
					LR	4181-0214-10	4121-0214-10
1st Bicuspid HK	-17	0	4	3	LL	4181-0224-11	4121-0224-11
					LR	4181-0214-11	4121-0214-11
2nd Bicuspid	-22	0	4	3	LL	4181-0225-10	4121-0225-10
					LR	4181-0215-10	4121-0215-10
2nd Bicuspid HK	-22	0	4	3	LL	4181-0225-11	4121-0225-11
					LR	4181-0215-11	4121-0215-11

MBT* Prescription

Maxillary	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Central	17	4	0	3.9	UL	4182-0121-10	4122-0121-10
					UR	4182-0111-10	4122-0111-10
Lateral	10	8	0	3.3	UL	4182-0122-10	4122-0122-10
					UR	4182-0112-10	4122-0112-10
Cuspid	-7	8	0	3.3	UL	4182-0123-00	4122-0123-00
					UR	4182-0113-00	4122-0113-00
Cuspid HK	-7	8	0	3.3	UR	4182-0113-01	4122-0113-01
					UL	4182-0123-01	4122-0123-01
Cuspid	0	8	0	3.3	UR	4182-0113-10	4122-0113-10
	0	8	0	3.3	UL	4182-0123-10	4122-0123-10
Cuspid HK	0	8	0	3.3	UL	4182-0123-11	4122-0123-11
	0	8	0	3.3	UR	4182-0113-11	4122-0113-11
1st & 2nd Bicuspid	-7	0	0	3.5	U	4182-0104-10	4122-0104-10
1st & 2nd Bicuspid HK	-7	0	0	3.5	UR	4182-0114-11	4122-0114-11
					UL	4182-0124-11	4122-0124-11
Kit						4182-000K-H*	4122-000K-H*

*indicate Hook on 3,45, 345 when ordering with hooks

MBT* Prescription

Mandibular	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Low Anterior	-6	0	0	2.6	U	4182-0201-10	4122-0201-10
Cuspid	-6	3	0	3.5	LR	4182-0213-00	4122-0213-00
					LL	4182-0223-00	4122-0223-00
Cuspid HK	-6	3	0	3.5	LR	4182-0213-01	4122-0213-01
					LL	4182-0223-01	4122-0223-01
Cuspid	0	3	0	3.5	LR	4182-0213-10	4122-0213-10
					LL	4182-0223-10	4122-0223-10
Cuspid HK	0	3	0	3.5	LR	4182-0213-11	4122-0213-11
					LL	4182-0223-11	4122-0223-11
1st Bicuspid	-12	0	0	3.5	U	4182-0204-10	4122-0204-10
1st Bicuspid HK	-12	0	0	3.5	LR	4182-0214-11	4122-0214-11
					LL	4182-0224-11	4122-0224-11
2nd Bicuspid	-17	0	0	3.5	U	4182-0205-10	4122-0205-10
2nd Bicuspid HK	-17	0	0	3.5	LR	4182-0215-11	4122-0215-11

CARVING
SELF LIGATING BRACKET

FEATURES

- Twin design
- No special instruments needed
- Low profile
- Easy opening clip
- Torque-in-base

