

BUCCAL TUBES

EASY FIT LOW PROFILE TUBES

Maxillary 1st & 2nd Molar

Torque	Offset	RX	M/D	L/R	.018 1stMolar	.022 1st Molar	0.018 2nd Molar	0.022 2nd Molar
-10	8	Roth	4.2	R	6181-0116-01	6121-0116-01	6181-0117-01	6121-0117-01
				L	6181-0126-01	6121-0126-01	6181-0127-01	6121-0127-01
-14	10	MBT	4.2	R	6182-0116-01	6122-0116-01	6182-0117-01	6122-0117-01
				L	6182-0126-01	6122-0126-01	6182-0127-01	6122-0127-01
0	0	STD	4.2	UR/LL	6180-0116-01	6120-0116-01	6180-0117-01	6120-0117-01
				UL/LR	6180-0126-01	6120-0126-01	6180-0127-01	6120-0127-01

Mandibular 1st & 2nd Molar

Torque	Offset	RX	M/D	L/R	.018 1stMolar	.022 1st Molar	0.018 2nd Molar	0.022 2nd Molar
-25	5	Roth	4.2	R	6181-0216-01	6121-0216-01	6181-0217-01	6121-0217-01
				L	6181-0226-01	6121-0226-01	6181-0227-01	6121-0227-01
-20	0	MBT	4.2	R	6182-0216-01	6122-0216-01	6182-0217-01	6122-0217-01
				L	6182-0226-01	6122-0226-01	6182-0227-01	6122-0227-01
0	0	STD	4.2	UL/LR	6180-0126-01	6120-0126-01	6180-0127-01	6120-0127-01
				UR/LL	6180-0116-01	6120-0116-01	6180-0117-01	6120-0117-01

FEATURES

- Low profile file to avoid occlusal interference
- Full funneled entrance helps guide wire into buccal tube
- Compound contour base anatomically conforms to the tooth
- 80 Gauge Foil Mesh

