

**ROTH\* Prescription**

Maxillary	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Central	12	5	0	3.9	UR	4381-0111-10	4321-0111-10
					UL	4381-0121-10	4321-0121-10
Lateral	8	9	0	3.3	UR	4381-0112-10	4321-0112-10
					UL	4381-0122-10	4321-0122-10
Cuspid	0	11	0	3.4	UR	4381-0113-10	4321-0113-10
					UL	4381-0123-10	4321-0123-10
Cuspid HK	0	11	0	3.4	UR	4381-0113-11	4321-0113-11
					UL	4381-0123-11	4321-0123-11
1st & 2nd Bicuspid	-7	0	0	3.4	U	4381-0104-10	4321-0104-10
1st & 2nd Bicuspid HK	-7	0	0	3.4	UR	4381-0114-11	4321-0114-11
					UL	4381-0124-11	4321-0124-11
Kit						4381-000K-H*	4321-000K-H*

\*indicate Hook on 3,45, 345 when ordering with hooks, and \*\*M if ordering with Mesh

**ROTH\* Prescription**

Mandibular	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Anterior	0	0	0	2.6	U	4381-0201-10	4321-0201-10
Cuspid	-11	5	0	3.4	LR	4381-0213-10	4321-0213-10
					LL	4381-0223-10	4321-0223-10
Cuspid HK	-11	5	0	3.4	LR	4381-0213-11	4321-0213-11
					LL	4381-0223-11	4321-0223-11
1st Bicuspid	-17	0	0	3.4	L	4381-0204-10	4321-0204-10
1st Bicuspid HK	-17	0	0	3.4	LR	4381-0214-11	4321-0214-11
					LL	4381-0224-11	4321-0224-11
2nd Bicuspid	-22	0	0	3.4	U	4381-0205-10	4321-0205-10
2nd Bicuspid HK	-22	0	0	3.4	LR	4381-0215-11	4321-0215-11
					LL	4381-0225-11	4321-0225-11

**MBT\* Prescription**

Maxillary	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Central	17	4	0	3.9	UR	4382-0111-20	4322-0111-20
					UL	4322-0121-20	4322-0121-20
Lateral	10	8	0	3.3	UR	4382-0112-20	4322-0112-20
					UL	4382-0122-20	4322-0122-20
Cuspid	-7	8	0	3.4	UR	4382-0113-00	4322-0113-00
					UL	4382-0123-00	4322-0123-00
Cuspid HK	-7	8	0	3.4	UR	4382-0113-01	4322-0113-01
					UL	4382-0123-01	4322-0123-01
Cuspid	0	8	0	3.4	UR	4382-0113-10	4322-0113-10
					UL	4382-0123-10	4322-0123-10
Cuspid HK	0	8	0	3.4	UR	4382-0113-11	4322-0113-11
					UL	4382-0123-11	4322-0123-11
1st & 2nd Bicuspid	-7	0	0	3.4	U	4382-0104-10	4322-0104-10
1st & 2nd Bicuspid HK					UR	4382-0114-11	4322-0114-11
					UL	4382-0124-11	4322-0124-11
Kit						4382-000K-H*	4322-000K-H*

\*indicate Hook on 3,45, 345 when ordering with hooks, and \*\*M if ordering with Mesh

**MBT\* Prescription**

Mandibular	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Anterior	-6	0	0	2.6	U	4382-0201-00	4322-0201-00
Cuspid	-6	3	0	3.4	LR	4382-0213-00	4322-0213-00
					LL	4382-0223-00	4322-0223-00
Cuspid HK	-6	3	0	3.4	LR	4382-0213-01	4322-0213-01
					LL	4382-0223-01	4322-0223-01
Cuspid	0	3	0	3.4	LR	4382-0213-10	4322-0213-10
					LL	4382-0223-10	4322-0223-10
Cuspid HK	0	3	0	3.4	LR	4382-0213-11	4322-0213-11
					LL	4382-0223-11	4322-0223-11
1st Bicuspid	-12	0	0	3.4	L	4382-0204-20	4322-0204-20
1st Bicuspid HK	-12	0	0	3.4	LR	4382-0214-21	4322-0214-21
					LL	4382-0224-21	4322-0224-21
2nd Bicuspid	-17	0	0	3.4	U	4382-0205-20	4322-0205-20
2nd Bicuspid HK	-17	0	0	3.4	LR	4382-0215-21	4322-0215-21
					LL	4382-0225-21	4322-0225-21

**LEAF**  
LOW PROFILE TWIN BRACKET

**FEATURES**

- High quality stainless steel (17-4)
- Full rotational control
- Compound Contour base
- Torque-in-base
- Also available with mesh base\*\*



\*\* When ordering with a mesh base please indicate M after the order number.

