

**ANTERIOR REGION**


Delivers lowest possible heat activated force to the teeth

**BICUSPID REGION**

Delivers progressively increased force to region

**POSTERIOR REGION**

Heavy force, concentrates its greatest force in this region



**FEATURES**

- Differentiated force value
- Broad size
- Up to 10 weeks without loss of corrective potential of the archwires
- 10 per pack

Round Size	Force Levels in Grams			Upper	Lower
	Anterior	Bicuspid	Posterior		
.014"	30	45	60	11-20-014-3F	12-20-014-3F
.016"	50	60	90	11-20-016-3F	12-20-016-3F
.018"	60	75	100	11-20-018-3F	12-20-018-3F

Round Size	Force Levels in Grams			Upper	Lower
	Anterior	Bicuspid	Posterior		
.016x.016"	50	70	100	11-20-1616-3F	12-20-1616-3F
.016x.022"	60	100	150	11-20-1622-3F	12-20-1622-3F
.017x.025"	90	160	250	11-20-1725-3F	12-20-1725-3F
.018x.018"	60	80	150	11-20-1818-3F	12-20-1818-3F
.018x.025"	100	170	260	11-20-1825-3F	12-20-1825-3F
.019"x.025"	120	200	280	11-20-1925-3F	12-20-1925-3F
.020"x.020"	90	150	20	11-20-2020-3F	12-20-2020-3F

