

ROTH* Prescription

Maxillary	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Central	12	5	0	3.5	UR	4581-0111-10	4521-0111-10
	12	5	0	3.5	UL	4581-0121-10	4521-0121-10
Central Hi-Torque	22	5	0	3.5	UR	4581-0111-20	4521-0111-20
	22	5	0	3.5	UL	4581-0121-20	4521-0121-20
Lateral	8	9	0	3	UR	4581-0112-10	4521-0112-10
	8	9	0	3	UL	4581-0122-10	4521-0122-10
Lateral Hi-Torque	14	9	0	3	UR	4581-0112-20	4521-0112-20
	14	9	0	3	UL	4581-0122-20	4521-0122-20
Cuspid	0	11	0	3.3	UR	4581-0113-10	4521-0113-10
	0	11	0	3.3	UL	4581-0123-10	4521-0123-10
Cuspid HK	0	11	0	3.3	UR	4581-0113-11	4521-0113-11
	0	11	0	3.3	UL	4581-0123-11	4521-0123-11
1st & 2nd Bicuspid	-7	0	0	3	U	4581-0104-10	4521-0104-10
1st & 2nd Bicuspid HK	-7	0	0	3	UR	4581-0114-11	4521-0114-11
	-7	0	0	3	UL	4581-0124-11	4521-0124-11
Kit						4581-000K-H*	4521-000K-H*

*indicate Hook on 3,45, 345 when ordering with hooks, and **M if ordering with Mesh

ROTH* Prescription

Mandibular	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Anterior	0	0	0	2.5	U	4581-0201-10	4521-0201-10
Cuspid	-11	5	0	3	LR	4581-0213-10	4521-0213-10
	-11	5	0	3	LL	4581-0223-10	4521-0223-10
Cuspid HK	-11	5	0	3	LR	4581-0213-11	4521-0213-11
	-11	5	0	3	LL	4581-0223-11	4521-0223-11
1st Bicuspid	-17	0	0	3	U	4581-0204-10	4521-0204-10
1st Bicuspid HK	-17	0	0	3	LR	4581-0214-11	4521-0214-11
	-17	0	0	3	LL	4581-0224-11	4521-0224-11
2nd Bicuspid	-22	0	0	3	U	4581-0205-10	4521-0205-10
2nd Bicuspid HK	-22	0	0	3	LR	4581-0215-11	4521-0215-11
	-22	0	0	3	LL	4581-0225-11	4521-0225-11

MBT* Prescription

Maxillary	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Central	17	4	0	3.9	UR	4582-0111-00	4522-0111-00
	17	4	0	3.9	UL	4582-0121-00	4522-0121-00
Lateral	10	8	0	3.3	UR	4582-0112-00	4522-0112-00
	10	8	0	3.3	UL	4582-0122-00	4522-0122-00
Cuspid	-7	8	0	3.3	UR	4582-0113-00	4522-0113-00
	-7	8	0	3.3	UL	4582-0123-00	4522-0123-00
Cuspid HK	-7	8	0	3.3	UR	4582-0113-01	4522-0113-01
	-7	8	0	3.3	UL	4582-0123-01	4522-0123-01
1st & 2nd Bicuspid	-7	0	0	3.5	U	4582-0104-00	4522-0104-00
1st & 2nd Bicuspid HK	-7	0	0	3.5	UR	4582-0114-01	4522-0114-01
	-7	0	0	3.5	UL	4582-0124-01	4522-0124-01
Kit						4582-000K-H*	4522-000K-H*

*indicate Hook on 3,45, 345 when ordering with hooks, and **M if ordering with Mesh

MBT* Prescription

Mandibular	Torque	Angulation	Rotation	M/D	L/R	.018	.022
Anterior	-6	0	0	2.6	U	4582-0201-00	4522-0201-00
Cuspid	-6	3	0	3.5	LR	4582-0213-00	4522-0213-00
	-6	3	0	3.5	LL	4582-0223-00	4522-0223-00
Cuspid HK	-6	3	0	3.5	LR	4582-0213-01	4522-0213-01
	-6	3	0	3.5	LL	4582-0223-01	4522-0223-01
1st Bicuspid	-12	0	0	3.5	U	4582-0204-20	4522-0204-20
1st Bicuspid HK	-12	0	0	3.5	LR	4582-0214-21	4522-0214-21
	-12	0	0	3.5	LL	4582-0224-21	4522-0224-21
2nd Bicuspid	-17	0	0	3.5	U	4582-0205-20	4522-0205-20
2nd Bicuspid HK	-17	0	0	3.5	LR	4582-0215-21	4522-0215-21
	-17	0	0	3.5	LL	4582-0225-21	4522-0225-21

FEATURES

- High quality stainless steel (17-4)
- One piece design to eliminate pad separation
- Low profile rhomboid design
- Torque-in-base
- Low friction bracket due to design
- Also available with mesh base**



** When ordering with a mesh base please indicate M after the order number.

