

ORDER INFORMATION

Lower	.018"	.020"	.018" x .025"	.019" x .025"
22 mm	12-10-018-22PS	12-10-020-22PS	12-10-1825-22PS	12-10-1925-22PS
24 mm	12-10-018-24PS	12-10-020-24PS	12-10-1825-24PS	12-10-1925-24PS
26 mm	12-10-018-26PS	12-10-020-26PS	12-10-1825-26PS	12-10-1925-26PS
28 mm	12-10-018-28PS	12-10-020-28PS	12-10-1825-28PS	12-10-1925-28PS
30 mm	12-10-018-30PS	12-10-020-30PS	12-10-1825-30PS	12-10-1925-30PS

Lower	.018"	.020"	.018" x .025"	.019" x .025"
32 mm	11-10-018-32PS	11-10-020-32PS	11-10-1825-32PS	11-10-1925-32PS
34 mm	11-10-018-34PS	11-10-020-34PS	11-10-1825-34PS	11-10-1925-34PS
36 mm	11-10-018-36PS	11-10-020-36PS	11-10-1825-36PS	11-10-1925-36PS
38 mm	11-10-018-38PS	11-10-020-38PS	11-10-1825-38PS	11-10-1925-38PS
40 mm	11-10-018-40PS	11-10-020-40PS	11-10-1825-40PS	11-10-1925-40PS



**STAINLESS STEEL
POSTED ARCHWIRES**

FEATURES TRUE AND EURO

- Great for retraction
- Solid brass post is easily formed
- Advanced brazing to eliminate separation of posts from archwire
- 10 per pack

